

# STERILE PETROLEUM JELLY

## SPECIFICATIONS

2 ml ampoule with 10 tips, amber glass,  
100 ampoules package

## COMPOSITION

100 % sterile petroleum jelly in conformity with European Pharmacopoeia

## USAGE

Sterile petroleum jelly is used to lubricate urinary probes particularly.

## PRECAUTIONS FOR USE

Not injectable.

To avoid any risk of contamination during placement of urinary probe, strict aseptic conditions should be accomplished.

Not to be used with latex probes. In fact, petroleum jelly alters latex porosity.

To keep out of reach of children.

Verify physical integrity of container to ensure sterile conditions prior to use.

## SHELF LIFE

5 years

<b>CDM LAVOISIER</b>		
<b>Laboratoires Chaix et Du Marais - 7, rue Labie - 75017 PARIS - FRANCE</b>		
<b>Tel : + 33 1 55 37 83 83</b>	<b>E-mail : <a href="mailto:contact@lavoisier.com">contact@lavoisier.com</a></b>	<b>Fax : + 33 1 55 37 83 84</b>